

Small, Simple, Smart ECG Patch

<sup>®</sup>  
**S-PATCH**



## Medical grade heart monitoring and analysis solution



Lightweight and convenient  
patch-type wearable device



Provide a comprehensive  
monitoring and analysis solution



# Innovative Lightweight & Comfortable Wearable

S-Patch is the Key  
to your seamless and unmatched  
continuous ECG monitoring.

## Holter

24/48 hour  
Holter

## CEM

Cardiac  
Event Monitor

## MCT

Mobile Cardiac  
Telemetry

## Extended Holter

14 days

**4-in-1 Test Capability**

All the remarkable power of **S-Patch**,

crafted into the lightness of just **two sheets of paper.**

## Recent S-Patch achievements

- S-Patch launched in over 16 countries
- S-Patch deployed to over 600 Hospitals and Clinics worldwide
- S-Patch received FDA 510K Clearance
- S-Patch certified CE, ISO 13485, MFDS

**GLOBAL  
CERTIFICATION**



# Single lead wearable patch

## Powered by Samsung Smart Health Processor



### Ultralight. All-day comfort.

- Weighing just 8g / 11g, designed to be forgotten



### Splashproof

- Water-resistant patch (IP55 grade)



### Cable-Free & Offline Test

- Built-in memory without data loss
- Offline mode with data storage
- Wireless data upload



### Reuse by quick swap

- Fast and simple replacement of battery & adhesives
- Compatible with commercially available coin battery and electrode

Clinical Accuracy.  
Continuous ECG.

## S-PATCH EX

Test Duration

**1~3 days**

Weight

**8g**

FDA clearance

**Aug. 2023**



## S-PATCH EXL

Test Duration

**14 days**

Weight

**11g**

FDA clearance

**Feb. 2024**


# Wear it. Be healthy.

## Transforming Healthcare from CURE to CARE

Wellysis is a digital healthcare solution company,  
a spin-off from Samsung SDS

S-Patch powered by Samsung Smart Health Processor was awarded  
HIMSS ASIAPAC INNOVATION AWARD FOR ARTIFICIAL INTELLIGENCE

If you want further information about Wellysis S-Patch solution,  
please contact [info@wellysis.com](mailto:info@wellysis.com), or access [www.wellysis.com](http://www.wellysis.com)



### Upcoming Features



Multi-Sensors



Personalization



Integration

[www.wellysis.com](http://www.wellysis.com)

© Wellysis | All rights reserved 2026

**Wellysis**